



# Pharmaceutical applications of UV/Visible Spectroscopy

Supervised by:

**Dr. Adel Magdy**

**Dr. Ahmed Hemdan**

Presented by team:

**Dr. Dev. Ibrahem Essam** : 120 11 345

**St. Mahmoud Hassan** : 11910463

**St. Mohamed Abd Elwahab** : 120 10 445

**St. Amr Hamdy** : 11910507

**St. Mohammed Abdelraheem Shabib** : 120 11 678



# ***Ultraviolet/Visible Spectroscopy ?***

- Ultraviolet/Visible (UV/Vis) light spectroscopy measures the absorbance of light by a sample in UV/Vis ranges to identify it.

## ***Advantages:***

- Easy to use: simple design of instrument.
- Fast results: quickly analyzes HPLC results [UV lamp].
- Maintains sample integrity: UV/Vis is a non-destructive technique compared to FTIR (Fourier-transform infrared spectroscopy).
- Highly sensitive in detecting organic compounds.

# ~: Applications :~

---

## ***1. Structure inspection of organic compounds:***

- Inspect organic molecules. From the location and combination of peaks, we find whether:
  - the compound is saturated or unsaturated.
  - hetero atoms are present or not....etc.

## ***2. Detection of Impurities:***

- Catalysts used in the reaction may still in the end product in small amounts as impurities.
- Additional peaks can be observed due to impurities in the sample compared with that of standard raw material.
- By also measuring the absorbance at  $\lambda_{\text{max}}$  of an impurity, the impurity can be detected.
- Ex: Benzene appears as a common impurity in cyclohexane. Its presence can be easily detected by its absorption at  $\lambda_{\text{max}}$  255 nm.

### ***3. Qualitative analysis:***

- Used for compounds that absorb UV/Vis radiation.
- Identification is done by comparing the absorption spectrum with the spectra of known compounds.
- Ex:
  - used for characterizing:
    - aromatic compounds.
    - aromatic olefins.

## 4. *Quantitative analysis:*

- Used for compounds that absorb UV//Vis radiation.
- This determination is based on Beer's law which is as follows:

$$A = \log I_0 / I_t = \log 1/T = -\log T = abc = \epsilon bc$$

Where:

$\epsilon$  is molar absorptivity.


$c$  is concentration.

$b$  is the length of the cuvette used in UV spectrophotometer.

- Other methods for quantitative analysis are as follows:
  - I. difference spectrophotometric method. [variable length]
  - II. calibration curve method. [used in labs]

## ***5. Dissociation constants of acids and bases:***

Consider the weak acid HB during its dissociation:



$$K_a = [\text{H}^+][\text{B}^-] / [\text{HB}]$$

$$-\log K_a = -\log ([\text{H}^+]) - \log ([\text{B}^-] / [\text{HB}])$$


$$\text{So, } \mathbf{pK_a = pH - \log ([B^-] / [HB])}$$

From the above,


$$\text{when } [\text{B}^-] = [\text{HB}] , \mathbf{pK_a = pH - \log ( 1 )}$$

and  $\log ( 1 ) = 0$  .so,  $\mathbf{pK_a = pH}$  at this point.

Measuring real-time changes in Absorbance of  $[HB]$  and  $[B^-]$  at  $\lambda_{\max}$  of each one of them.



The **pH** at this point is **pK<sub>a</sub>**





## **6. *Chemical kinetics:***

The UV radiation is passed through the reaction cell and the absorbance changes can be observed.

## **7. *As HPLC detector:***

- A UV/Vis spectrophotometer may be used as a detector for HPLC (high-performance liquid chromatography) (ie: UV lamp).
- The presence of an analyte gives a response which can be assumed to be proportional to the concentration.

## ***8. Molecular weight determination:***

- Molecular weights can be measured spectrophotometrically by preparing a suitable derivative of unknown compound.
- For example,
  - if we want to determine the molecular weight of amine:
    - it is converted into amine picrate.
    - known weight of amine picrate is dissolved in a litre of solution.
    - its optical density is measured at  $\lambda_{\text{max}}$  380 nm.

After this, the concentration of the solution in moles per litre (Molar) can be calculated by using the following formula:

- "c" can be calculated using above equation, the weight "w" of amine picrate is known.

$$C = \frac{\log I_0 / I_t}{\epsilon_{\max} \times l}$$

- From "c" and "w",  
molecular weight of amine picrate can be calculated:  
c = number of moles / volume in Liter  
"c" is known, volume is 1 Liter. So, moles can be calculated.  
And "w" represents these moles.  
So, mol.wt of amine picrate  
can be calculated by cross multiplication.
- And the molecular weight of amine can be calculated using the molecular weight of amine picrate.

## *~: References :~*

- Dr. Adel's handout. [;]
- Dr. Ahmed's inspiration.
- <https://www.pharmatutor.org/pharma-analysis/analytical-aspects-of-uv-visible-spectroscopy/applications.html>
- [https://www.researchgate.net/publication/322935261\\_Ultraviolet\\_spectroscopy\\_and\\_its\\_pharmaceutical\\_applications-\\_A\\_brief\\_review](https://www.researchgate.net/publication/322935261_Ultraviolet_spectroscopy_and_its_pharmaceutical_applications-_A_brief_review)
- <https://www.laboratory-equipment.com/blog/absorption-spectroscopy-in-pharmaceutical-analysis/>

*Thank  
You!*